

Raz Enterprises

INDEX

1. BIPOLAR FORCEPS	05-07
2. MONOPOLAR FORCEPS	08
3. NON-STICK BIPOLAR FORCEPS	09-10
4. BIPOLAR HAND-SWITCH FORCEPS	11
5. IRRIGATION BIPOLAR FORCEPS	12
6. MONOPOLAR HAND-SWITCH FORCEPS	13
7. EUROPEAN BIPOLAR FORCEPS	14-16
8. EUROPEAN NON-STICK BIPOLAR FORCEPS	17-18
9. EUROPEAN IRRIGATION BIPOLAR FORCEPS	19
10. DISPOSABLE BIPOLAR FORCEPS	20-21
11. DIATHERMY INSTRUMENTS	22-23
12. BIPOLAR ARTERY SEALER	24-25
13. BIPOLAR ELECTRODES	26
14. BIPOLAR ELECTRODES (SINGLE USE)	27
15. BIPOLAR CABLE	29-30
16. MONOPOLAR CABLE	31
17. ELECTROSURGICAL PENCIL	33
18. ELECTRODES (2.4 mm)	34
19. LOOP ELECTRODES (2.4 mm)	35
20. ELECTRODES (4.0 mm)	36
21. LOOP ELECTRODES (4.0 mm)	37
22. PATIENT RETURN PLATE	38
23. PATIENT RETURN PLATE (SILICONE)	39
24. LAPAROSCOPIC ELECTRODES	40
25. ARTHROSCOPIC ELECTRODES	41
26. ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY	43-46
27. LASER SURGERY INSTRUMENTS FOR GYNECOLOGY	47-48





OPHTHALMIC BIPOLAR FORCEPS



McPHERSON
(STR.)
110-100



McPHERSON
(CRV.)
110-110



COAPTATION
NADLER
110-120



110-130



IRIS (STR.)
110-180



IRIS (CRV.)
110-190



JEWELER #5
110-140



JEWELER #4
110-150



TENZEL
110-160



JEWELER #7
110-170

Cat #	Description	Tip (mm)	Size
110-100	McPHERSON (STR.)	0.5	3 1/2" (9.0 cm)
110-110	McPHERSON (CRV.)	0.5	3 1/2" (9.0 cm)
110-120	COAPTATION NADLER	0.5	4" (10.2 cm)
110-130	COAPTATION	0.5	4 1/4" (10.8 cm)
110-140	JEWELER #5	MICRO, 0.5	4 1/2" (11.5 cm)
110-150	JEWELER #4	MICRO, 0.5	4 1/2" (11.5 cm)
110-160	TENZEL	MICRO, 0.5	4" (10.2 cm)
110-170	JEWELER #7	MICRO, 0.5	4 1/2" (11.5 cm)
110-180	IRIS (STR.)	1.0	4 1/4" (10.8 cm)
110-190	IRIS (CRV.)	1.0	4 1/4" (10.8 cm)

BIPOLAR FORCEPS



NEEDLE 0.5mm 0.7mm 1.0mm 1.5mm 2.0mm
Available Tip Widths Always Specify Tips Width

For your convenience, add the following tip extension to product code.
Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) and 2.0 mm (-20)



PLASTIC & RECONSTRUCTIVE

BIPOLAR FORCEPS

ADSON
110-200



SEMKIN
110-220



CUSHING
110-230



CUSHING
110-235



CUSHING (GURSH)
110-240



GERALD (STR.)
110-300



GERALD (CRV.)
110-310



GERALD BAYONET
110-320



For your convenience, add the following tip extension to product code.

Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) and 2.0 mm (-20)



NON-STICK BIPOLAR FORCEPS

No more undesirable tissue sticking during coagulation procedure. Thanks to RAZ new **Ag-Tech** non-stick technology, the forceps stay clean- virtually eliminating the need for cleaning during coagulation. By applying a material with excellent thermal properties for the tips which prevents the sticking of tissues, the non stick effect is permanent and will not reduce with the use of forceps.

MCPHERSON (STR.)
110-100 NS



MCPHERSON (CRV.)
110-110 NS



JEWELER #5
110-140 NS



JEWELER #4
110-150 NS



ADSON
110-200 NS



SEMKN
110-220 NS



CUSHING
110-230 NS



CUSHING
110-235 NS





NON-STICK BIPOLAR FORCEPS



GERALD BAYONET
110-320 NS



JANSEN BAYONET
110-360 NS



SCOVILLE-GREENWOOD
110-430 NS



YASARGIL
110-440 NS



HARDY BAYONET
110-460 NS



BIPOLAR HAND-SWITCH FORCEPS


JEWELER #4
110-150 HS




SEMKIN
110-220 HS



JANSEN BAYONET
110-370 HS



SCOVILLE-GREENWOOD
110-430 HS



More Bipolar Hand-switch forceps are available.
Use "HS" extension to any regular bipolar forceps.



IRRIGATION BIPOLAR FORCEPS



HARDY BAYONET
110-460 IR



JANSEN BAYONET
110-370 IR



SCOVILLE-GREENWOOD
110-430 IR



GERALD BAYONET
110-320 IR



CUSHING
110-230 IR

More Bipolar Irrigation forceps are available.
Use "IR" extension to any regular bipolar forceps.



MONOPOLAR HAND-SWITCH FORCEPS



JANSEN BAYONET
110-560 HS



HAND ACTUATED
SWITCH

JANSEN BAYONET
110-580 HS



POTTS-SMITH
110-540 HS



POTTS-SMITH
110-550 HS





EUROPEAN BIPOLAR FORCEPS

Raz Enterprises

P.O. BOX NO. 449.

SIALKOT-51310-Pakistan.

Phone: +92-52-3551857

Fax: +92-52-3558041

E-mail: Info@Raz-Enterprises.com



MCPHERSON
310-100



MCPHERSON (ANGLED)
310-102



JEWELER #5
310-110



JEWELER #5 (ANGLED)
310-112



JEWELER #7
310-115



ADSON
310-300



NEEDLE 0.5mm 0.7mm 1.0mm 1.5mm 2.0mm
Available Tip Widths Always Specify Tips Width

For your convenience, add the following tip extension to product code.

Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) and 2.0 mm (-20)

HIGH FREQUENCY INNOVATION



EUROPEAN BIPOLAR FORCEPS

STRAIGHT
310-120



STRAIGHT
310-130



STRAIGHT
310-140



STRAIGHT
310-150



STRAIGHT (ANGLED)
310-122



STRAIGHT (ANGLED)
310-132



STRAIGHT (ANGLED)
310-142



STRAIGHT (ANGLED)
310-152



For your convenience, add the following tip extension to product code.

Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) and 2.0 mm (-20)



EUROPEAN BIPOLAR FORCEPS

Cat #	Description	Tip (mm)	Size
310-160	BAYONET	0.5, 1.0, 2.0	6 1/2" (16.5 cm)
310-170	BAYONET	0.5, 1.0, 2.0	7" (18.0 cm)
310-180	BAYONET	0.5, 1.0, 2.0	7 3/4" (20.0 cm)
310-190	BAYONET	0.5, 1.0, 2.0	8 3/4" (22.0 cm)
310-200	BAYONET	0.5, 1.0, 2.0	9 1/2" (24.0 cm)
310-162	BAYONET (ANGLED UP)	0.5, 1.0, 2.0	6 1/2" (16.5 cm)
310-172	BAYONET (ANGLED UP)	0.5, 1.0, 2.0	7" (18.0 cm)
310-182	BAYONET (ANGLED UP)	0.5, 1.0, 2.0	7 3/4" (20.0 cm)

BAYONET
310-160



BAYONET
310-170



BAYONET
310-180



BAYONET (ANGLED UP)
310-162



BAYONET (ANGLED UP)
310-172



BAYONET (ANGLED UP)
310-182



BAYONET (ANGLED DOWN)
310-165



BAYONET (ANGLED DOWN)
310-175



BAYONET (ANGLED DOWN)
310-185



For your convenience, add the following tip extension to product code.

Needle Tip (-00), 0.5 mm (-05), 1.0 mm (-10) and 2.0 mm (-20)



EUROPEAN NON-STICK BIPOLAR FORCEPS



McPHERSON
310-100 NS



McPHERSON (ANGLED)
310-102 NS



JEWELER #5
310-110 NS



JEWELER #5 (ANGLED)
310-112 NS

No more undesirable tissue sticking during coagulation procedure. Thanks to RAZ new **Ag-Tech** non-stick technology, the forceps stay clean- virtually eliminating the need for cleaning during coagulation. By applying a material with excellent thermal properties for the tips which prevents the sticking of tissues, the non stick effect is permanent and will not reduce with the use of forceps.



ADSON
310-300 NS



EUROPEAN NON-STICK BIPOLAR FORCEPS

STRAIGHT
310-120 NS



STRAIGHT
310-130 NS



STRAIGHT
310-140 NS



STRAIGHT
310-150 NS



BAYONET
310-160 NS



BAYONET
310-170 NS



BAYONET
310-180 NS



BAYONET
310-190 NS





EUROPEAN IRRIGATION


BIPOLAR FORCEPS



BAYONET
310-160-IR



BAYONET
310-180-IR



BAYONET
310-150-IR



BAYONET
310-120-IR



DISPOSABLE BIPOLAR FORCEPS



GERALD (Str)
210-300



JANSEN BAYONET
(IRRIGATION)
210-360 IR



JANSEN BAYONET
210-370



JEWELER #5
210-140



ADSON
210-200



DISPOSABLE BIPOLAR CABLE
190-100 D





DISPOSABLE BIPOLAR FORCEPS



DISPOSABLE BIPOLAR CABLE



JEWELER #5
410-110



STRAIGHT FORCEPS
410-140



BAYONET FORCEPS
410-180



STRAIGHT FORCEPS
(Disposable without cable)
510-140



STRAIGHT ANG
(Disposable without cable)
510-142



DIATHERMY INSTRUMENTS



TUNER WARWICK (STR.)
320-200



TUNER WARWICK (CRV.)
320-202



SURGICAL SCISSORS (STR.)
320-220



SURGICAL SCISSORS (CRV.)
320-222



WILSON HEY FORCEPS (STR.)
320-240



WILSON HEY FORCEPS (CRV.)
320-242



ARTERY FORCEPS (STR.)
320-260



ARTERY FORCEPS (CRV.)
320-262



ABBEY NEEDLE (17.5 cm)
320-190



RICHIE FORCEPS (Str.)
320-290

ENGLISH MONOPOLAR
(MALE 4.8 mm CONNECTOR)



EUROPEAN MONOPOLAR
(FEMALE 4.0 mm CONNECTOR)



Cat #	Description	Size	Cat #	Description	Size
320-200	TUNER-WARWICK SCISSORS (STR.)	7" (17.8 cm)	330-200	TUNER-WARWICK SCISSORS (STR.)	7" (17.8 cm)
320-202	TUNER-WARWICK SCISSORS (CRV.)	7" (17.8 cm)	330-202	TUNER-WARWICK SCISSORS (CRV.)	7" (17.8 cm)
320-210	TUNER-WARWICK SCISSORS (STR.)	9" (23.0 cm)	330-210	TUNER-WARWICK SCISSORS (STR.)	9" (23.0 cm)
320-212	TUNER-WARWICK SCISSORS (CRV.)	9" (23.0 cm)	330-212	TUNER-WARWICK SCISSORS (CRV.)	9" (23.0 cm)
320-220	SURGICAL SCISSORS (STR.)	7" (17.8 cm)	330-220	SURGICAL SCISSORS (STR.)	7" (17.8 cm)
320-222	SURGICAL SCISSORS (CRV.)	7" (17.8 cm)	330-222	SURGICAL SCISSORS (CRV.)	7" (17.8 cm)
320-230	SURGICAL SCISSORS (STR.)	9" (23.0 cm)	330-230	SURGICAL SCISSORS (STR.)	9" (23.0 cm)
320-232	SURGICAL SCISSORS (CRV.)	9" (23.0 cm)	330-232	SURGICAL SCISSORS (CRV.)	9" (23.0 cm)
320-240	WILSON HEY FORCEPS (STR.)	7" (17.8 cm)	330-240	WILSON HEY FORCEPS (STR.)	7" (17.8 cm)
320-242	WILSON HEY FORCEPS (CRV.)	7" (17.8 cm)	330-242	WILSON HEY FORCEPS (CRV.)	7" (17.8 cm)
320-250	WILSON HEY FORCEPS (STR.)	9" (23.0 cm)	330-250	WILSON HEY FORCEPS (STR.)	9" (23.0 cm)
320-252	WILSON HEY FORCEPS (CRV.)	9" (23.0 cm)	330-252	WILSON HEY FORCEPS (CRV.)	9" (23.0 cm)
320-260	ARTERY FORCEPS (STR.)	7" (17.8 cm)	330-260	ARTERY FORCEPS (STR.)	7" (17.8 cm)
320-262	ARTERY FORCEPS (CRV.)	7" (17.8 cm)	330-262	ARTERY FORCEPS (CRV.)	7" (17.8 cm)
320-270	ARTERY FORCEPS (STR.)	9" (23.0 cm)	330-270	ARTERY FORCEPS (STR.)	9" (23.0 cm)
320-272	ARTERY FORCEPS (CRV.)	9" (23.0 cm)	330-272	ARTERY FORCEPS (CRV.)	9" (23.0 cm)
320-290	RICHIE FORCEPS	7" (17.8 cm)	330-290	RICHIE FORCEPS	7" (17.8 cm)
320-292	RICHIE FORCEPS	8" (20.0 cm)	330-292	RICHIE FORCEPS	8" (20.0 cm)
320-296	RICHIE FORCEPS	9" (23.0 cm)	330-296	RICHIE FORCEPS	9" (23.0 cm)



DIATHERMY INSTRUMENTS



ADSON
320-100



McINDOE
320-110



GILLIES
320-120



DISSECTING
320-130



GERALD (1x2)
320-140



ADSON BAYONET
320-150



ADSON BAYONET (ANG. UP)
320-160



REMINGTON HOBBS
320-170

ENGLISH MONOPOLAR (MALE 4.8 mm CONNECTOR)

4.8 mm



Cat #	Description	Tip (mm)	Size
320-100	ADSON	1.0	4 ^{3/4} " (12.0 cm)
320-110	McINDOE	1.0, 2.0	6" (15.0 cm)
320-120	GILLIES (1X2)	—	6" (15.0 cm)
320-122	GILLIES (1X2)	—	7" (17.8 cm)
320-124	GILLIES (1X2)	—	8" (20.0 cm)
320-126	GILLIES (1X2)	—	9" (24.0 cm)
320-130	DISSECTING FORCEPS	1.0, 2.0	6" (15.0 cm)
320-132	DISSECTING FORCEPS	1.0, 2.0	7" (17.8 cm)
320-134	DISSECTING FORCEPS	1.0, 2.0	8 ^{1/4} " (21.0 cm)
320-136	DISSECTING FORCEPS	1.0, 2.0	9 ^{1/2} " (24.1 cm)
320-140	GERALD (1X2)	—	7" (17.8 cm)
320-150	ADSON BAYONET	1.0, 2.0	6 ^{1/4} " (16.0 cm)
320-152	ADSON BAYONET	1.0, 2.0	7 ^{1/2} " (19.1 cm)
320-160	ADSON BAYONET (ANG. UP)	1.0, 2.0	6 ^{1/4} " (16.0 cm)
320-162	ADSON BAYONET (ANG. UP)	1.0, 2.0	7 ^{1/2} " (19.1 cm)
320-164	ADSON BAYONET (ANG. UP)	1.0, 2.0	8 ^{1/4} " (21.0 cm)
320-170	REMINGTON HOB	1.0, 2.0	9 ^{1/2} " (24.1 cm)

EUROPEAN MONOPOLAR (FEMALE 4.0 mm CONNECTOR)

4.0 mm



Cat #	Description	Tip (mm)	Size
330-100	ADSON	1.0	4 ^{3/4} " (12.0 cm)
330-110	McINDOE	1.0, 2.0	6" (15.0 cm)
330-120	GILLIES (1X2)	—	6" (15.0 cm)
330-122	GILLIES (1X2)	—	7" (17.8 cm)
330-124	GILLIES (1X2)	—	8" (20.0 cm)
330-126	GILLIES (1X2)	—	9" (24.0 cm)
330-130	DISSECTING FORCEPS	1.0, 2.0	6" (15.0 cm)
330-132	DISSECTING FORCEPS	1.0, 2.0	7" (17.8 cm)
330-134	DISSECTING FORCEPS	1.0, 2.0	8 ^{1/4} " (21.0 cm)
330-136	DISSECTING FORCEPS	1.0, 2.0	9 ^{1/2} " (24.1 cm)
330-140	GERALD (1X2)	—	7" (17.8 cm)
330-150	ADSON BAYONET	1.0, 2.0	6 ^{1/4} " (16.0 cm)
330-152	ADSON BAYONET	1.0, 2.0	7 ^{1/2} " (19.1 cm)
330-160	ADSON BAYONET (ANG. UP)	1.0, 2.0	6 ^{1/4} " (16.0 cm)
330-162	ADSON BAYONET (ANG. UP)	1.0, 2.0	7 ^{1/2} " (19.1 cm)
330-164	ADSON BAYONET (ANG. UP)	1.0, 2.0	8 ^{1/4} " (21.0 cm)
330-170	REMINGTON HOB	1.0, 2.0	9 ^{1/2} " (24.1 cm)

EUROPEAN MONOPOLAR (MALE 4.0 mm CONNECTOR)

4.0 mm



Cat #	Description	Tip (mm)	Size
340-100	ADSON	1.0	4 ^{3/4} " (12.0 cm)
340-110	McINDOE	1.0, 2.0	6" (15.0 cm)
340-120	GILLIES (1X2)	—	6" (15.0 cm)
340-122	GILLIES (1X2)	—	7" (17.8 cm)
340-124	GILLIES (1X2)	—	8" (20.0 cm)
340-126	GILLIES (1X2)	—	9" (24.0 cm)
340-130	DISSECTING FORCEPS	1.0, 2.0	6" (15.0 cm)
340-132	DISSECTING FORCEPS	1.0, 2.0	7" (17.8 cm)
340-134	DISSECTING FORCEPS	1.0, 2.0	8 ^{1/4} " (21.0 cm)
340-136	DISSECTING FORCEPS	1.0, 2.0	9 ^{1/2} " (24.1 cm)
340-140	GERALD (1X2)	—	7" (17.8 cm)
340-150	ADSON BAYONET	1.0, 2.0	6 ^{1/4} " (16.0 cm)
340-152	ADSON BAYONET	1.0, 2.0	7 ^{1/2} " (19.1 cm)
340-160	ADSON BAYONET (ANG. UP)	1.0, 2.0	6 ^{1/4} " (16.0 cm)
340-162	ADSON BAYONET (ANG. UP)	1.0, 2.0	7 ^{1/2} " (19.1 cm)
340-164	ADSON BAYONET (ANG. UP)	1.0, 2.0	8 ^{1/4} " (21.0 cm)
340-170	REMINGTON HOB	1.0, 2.0	9 ^{1/2} " (24.1 cm)

MONOPOLAR (WITHOUT CONNECTOR)




Cat #	Description	Tip (mm)	Size
350-100	ADSON	1.0	4 ^{3/4} " (12.0 cm)
350-110	McINDOE	1.0, 2.0	6" (15.0 cm)
350-120	GILLIES (1X2)	—	6" (15.0 cm)
350-122	GILLIES (1X2)	—	7" (17.8 cm)
350-124	GILLIES (1X2)	—	8" (20.0 cm)
350-126	GILLIES (1X2)	—	9" (24.0 cm)
350-130	DISSECTING FORCEPS	1.0, 2.0	6" (15.0 cm)
350-132	DISSECTING FORCEPS	1.0, 2.0	7" (17.8 cm)
350-134	DISSECTING FORCEPS	1.0, 2.0	8 ^{1/4} " (21.0 cm)
350-136	DISSECTING FORCEPS	1.0, 2.0	9 ^{1/2} " (24.1 cm)
350-140	GERALD (1X2)	—	7" (17.8 cm)
350-150	ADSON BAYONET	1.0, 2.0	6 ^{1/4} " (16.0 cm)
350-152	ADSON BAYONET	1.0, 2.0	7 ^{1/2} " (19.1 cm)
350-160	ADSON BAYONET (ANG. UP)	1.0, 2.0	6 ^{1/4} " (16.0 cm)
350-162	ADSON BAYONET (ANG. UP)	1.0, 2.0	7 ^{1/2} " (19.1 cm)
350-164	ADSON BAYONET (ANG. UP)	1.0, 2.0	8 ^{1/4} " (21.0 cm)
350-170	REMINGTON HOB	1.0, 2.0	9 ^{1/2} " (24.1 cm)

HIGH FREQUENCY INNOVATION



BIPOLAR ARTERY SEALER



DETACHABLE CONNECTING CABLE
190-180

HIGH FREQUENCY INNOVATION



BIPOLAR ARTERY SEALER



ARTERY SEALER
110-760



ARTERY SEALER
110-764



BIPOLAR ELECTRODES

Cat #	Description	Size
310-500	STRAIGHT	18.0 cm
310-510	STRAIGHT	20.0 cm
310-520	STRAIGHT	22.0 cm
310-530	STRAIGHT	27.0 cm
310-540	ANGLED	18.0 cm
310-550	ANGLED	20.0 cm
310-560	ANGLED	22.0 cm
310-570	ANGLED (UPWARD TIP)	16.5 cm
310-580	ANGLED (UPWARD TIP)	20.0 cm
310-590	ANGLED 45°	19.0 cm
310-600	ANGLED 90°	16.0 cm



STRAIGHT
310-500



ANGLED
(Upward Tip)
310-570



ANGLED
310-540




ANGLED 45°
310-590



ANGLED 90°
310-600



BI-TURBINATE PROBE (17.5 cm)
110-700



HIGH FREQUENCY INNOVATION



BIPOLAR ELECTRODES

SINGLE USE

Cat #	Description	Size
510-500	STRAIGHT	18.0 cm
510-510	STRAIGHT	20.0 cm
510-520	STRAIGHT	22.0 cm
510-530	STRAIGHT	27.0 cm
510-540	ANGLED	18.0 cm
510-550	ANGLED	20.0 cm
510-560	ANGLED	22.0 cm
510-570	ANGLED (UPWARD TIP)	16.5 cm
510-580	ANGLED (UPWARD TIP)	20.0 cm
510-590	ANGLED 45°	19.0 cm
510-600	ANGLED 90°	16.0 cm



STRAIGHT
510-500

ANGLED
(Upward Tip)
510-570

ANGLED
510-540

ANGLED 45°
510-590

ANGLED 90°
510-600





BIPOLAR CABLES



SILICONE BIPOLAR CABLE 10' (3.0 M)

190-100

Instrument Connector: USA 2-Pin Plug

Generator Connector: 4mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus
& Most Bipolar Generator



SILICONE EUROPEAN BIPOLAR CABLE 10' (3.0 M)

190-110

Instrument Connector: European Flat Plug

Generator Connector: 4mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus
& Most Bipolar Generator



SILICONE EUROPEAN BIPOLAR CABLE 10' (3.0 M)

190-120

Instrument Connector: European Flat Plug

Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Berchtold, Wolf




SILICONE BIPOLAR CABLE 10' (3.0 M)

190-130

Instrument Connector: USA 2-Pin Plug

Generator Connector: Martin Type

Compatibility: Aesculap, Martin, Berchtold, Wolf



SILICONE EUROPEAN BIPOLAR CABLE 10' (3.0 M)

190-140

Instrument Connector: European Flat Plug


Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens

CABLES



BIPOLAR CABLES




SILICONE BIPOLAR CABLE 10' (3.0 M)

190-150

Instrument Connector: USA 2-Pin Plug

Generator Connector: Erbe Type

Compatibility: Erbe, Downs, Siemens




SILICONE EUROPEAN BIPOLAR CABLE 10' (3.0 M)

190-160

Instrument Connector: European Flat Plug

Generator Connector: Pin Distance 29mm

Compatibility: Valleylab, EMC, Dolley



SILICONE BIPOLAR CABLE 10' (3.0 M)

190-170

Instrument Connector: USA 2-Pin Plug

Generator Connector: Pin Distance 29mm

Compatibility: Valleylab, EMC, Dolley



SILICONE BIPOLAR CABLE 10' (3.0 M)

190-180

Instrument Connector: Bipolar Artery Sealer

Generator Connector: 4mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus
& Most Bipolar Generator



SILICONE BIPOLAR CABLE 10' (3.0 M)

190-185

Instrument Connector: Ethicon

Generator Connector: 4mm Banana Plug

Compatibility: Valleylab, Conmed, Codman, Olympus
& Most Bipolar Generator



MONOPOLAR CABLES



SILICONE MONOPOLAR CABLE 10' (3.0 M)

190-230

Instrument Connector: Screw Type

Generator Connector: Standard Bovie® 8mm Plug

Compatibility: Valleylab, Conmed, Bovie®



SILICONE MONOPOLAR ENDOSCOPIC CABLE 10' (3.0 M)

190-240

Instrument Connector: 4mm Female Fitting

Generator Connector: Standard Bovie® 8mm Plug

Compatibility: Valleylab, Conmed, Bovie®



SILICONE MONOPOLAR ENDOSCOPIC CABLE 10' (3.0 M)

190-250

Instrument Connector: 4mm Female Fitting

Generator Connector: Standard 4mm Banana Plug

Compatibility: Most ES Generators



SILICONE DIATHERMY CABLE 10' (3.0 M)

190-260

Instrument Connector: 4mm Male Fitting (European Diathermy)

Generator Connector: Standard 4mm Banana Plug

Compatibility: Most ES Generators



SILICONE DIATHERMY CABLE 10' (3.0 M)

190-280

Instrument Connector:

4.8mm Female Fitting(English Diathermy)

Generator Connector:

Standard 4mm Banana Plug

Compatibility: Most ES Generators



SILICONE DIATHERMY CABLE 10' (3.0 M)

190-270

Instrument Connector:

4.8mm Female Fitting(English Diathermy)

Generator Connector:

Standard Bovie® 8mm Plug

Compatibility: Valleylab, Conmed, Bovie®





ROCKER-SWITCH PENCIL

SINGLE USE

150-104 (3m Cable)



FOOT CONTROL PENCIL

SINGLE USE

150-102 (3m Cable)

150-103 (5m Cable)



FINGER-SWITCH PENCIL

REUSABLE

150-105 (3m Silicon Cable)



ELECTRODE HANDLE

REUSABLE

150-107 (3m Silicon Cable)



FINGER-SWITCH PENCIL

SINGLE USE

150-100 (3m Cable)

150-101 (5m Cable)

ELECTROSURGICAL PENCIL



ELECTRODES (2.4mm)



BLADE ELECTRODE

Disposable: 150-110, Reusable: 152-110



BLADE ELECTRODE (CRV.)

Disposable: 150-112, Reusable: 152-112



NEEDLE ELECTRODE

Disposable: 150-120, Reusable: 152-120



NEEDLE ELECTRODE (CRV.)

Disposable: 150-122, Reusable: 152-122



BALL ELECTRODE (2mm)

Disposable: 150-130, Reusable: 152-130



BALL ELECTRODE CRV. (2mm)

Disposable: 150-132, Reusable: 152-132



BALL ELECTRODE (3mm)

Disposable: 150-140, Reusable: 152-140



BALL ELECTRODE CRV. (3mm)

Disposable: 150-142, Reusable: 152-142



BALL ELECTRODE (4mm)

Disposable: 150-150, Reusable: 152-150



BALL ELECTRODE CRV. (4mm)

Disposable: 150-152, Reusable: 152-152



BALL ELECTRODE (5mm)

Disposable: 150-160, Reusable: 152-160



BALL ELECTRODE CRV. (5mm)

Disposable: 150-162, Reusable: 152-162



BLADE ELECTRODE (Shaft 13.0cm)

Disposable: 150-115, Reusable: 152-115



NEEDLE ELECTRODE (Shaft 13.0cm)

Disposable: 150-125, Reusable: 152-125




BALL ELECTRODE 3mm (Shaft 13.0cm)

Disposable: 150-145, Reusable: 152-145



BALL ELECTRODE 5mm (Shaft 13.0cm)

Disposable: 150-165, Reusable: 152-165



CONIZATION ELECTRODE

Disposable: 150-195, Reusable: 152-195



CIRCLE LOOP ELECTRODE (4.0mm)

Reusable: 152-510



CIRCLE LOOP ELECTRODE (6.0mm)

Reusable: 152-520



CIRCLE LOOP ELECTRODE (8.0mm)

Reusable: 152-530



DIAMOND ELECTRODE

Reusable: 152-540



FINE NEEDLE ELECTRODE

Reusable: 152-550



TIP CLEANER 150-109

**Raz
Enterprises**

P.O. BOX NO. 449.
SIALKOT-51310-Pakistan.
Phone: +92-52-3551857
Fax: +92-52-3558041
E-mail: Info@Raz-Enterprises.com



HOLSTER 150-108

HIGH FREQUENCY INNOVATION



LOOP ELECTRODES (2.4mm)

Raz Enterprises

P.O. BOX NO. 449.

SIALKOT-51310-Pakistan.

Phone: +92-52-3551857

Fax: +92-52-3558041

E-mail: Info@Raz-Enterprises.com



Size	1.0X0.4 cm	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm
Disposable	150-170-04	150-170-05	150-170-07	150-170-08	150-170-10
Reusable	152-170-04	152-170-05	152-170-07	152-170-08	152-170-10



Size	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm	1.0X1.5 cm
Disposable	150-175-05	150-175-07	150-175-08	150-175-10	150-175-15
Reusable	152-175-05	152-175-07	152-175-08	152-175-10	152-175-15



Size	1.5X0.5 cm	1.5X0.7 cm	1.5X0.8 cm	1.5X1.0 cm	1.5X1.2 cm	1.5X1.5 cm
Disposable	150-180-05	150-180-07	150-180-08	150-180-10	150-180-12	150-180-15
Reusable	152-180-05	152-180-07	152-180-08	152-180-10	152-180-12	152-180-15



Size	2.0X0.8 cm	2.0X1.0 cm	2.0X1.2 cm	2.0X1.3 cm	2.0X1.5 cm	2.0X2.0 cm
Disposable	150-190-08	150-190-10	150-190-12	150-190-13	150-190-15	150-190-20
Reusable	152-190-08	152-190-10	152-190-12	152-190-13	152-190-15	152-190-20

HIGH FREQUENCY INNOVATION



ELECTRODES (4.0mm)



BLADE ELECTRODE
Reusable 162-110



BLADE ELECTRODE (CRV.)
Reusable 162-112



NEEDLE ELECTRODE
Reusable 162-120



NEEDLE ELECTRODE (CRV.)
Reusable 162-122



BALL ELECTRODE (2mm)
Reusable 162-130



BALL ELECTRODE (CRV.) (2mm)
Reusable 162-132



BALL ELECTRODE (3mm)
Reusable 162-140



BALL ELECTRODE (CRV.) (3mm)
Reusable 162-142



BALL ELECTRODE (4mm)
Reusable 162-150



BALL ELECTRODE (CRV.) (4mm)
Reusable 162-152



BALL ELECTRODE (5mm)
Reusable 162-160



BALL ELECTRODE (CRV.) (5mm)
Reusable 162-162



BLADE ELECTRODE (Shaft 13.0 cm)
Reusable 162-115



NEEDLE ELECTRODE (Shaft 13.0 cm)
Reusable 162-125



BALL ELECTRODE 3mm (Shaft 13.0 cm)
Reusable 162-145



BALL ELECTRODE 4mm (Shaft 13.0 cm)
Reusable 162-155



BALL ELECTRODE 5mm (Shaft 13.0 cm)
Reusable 162-165



ROUND LOOP ELECTRODE (7.0cm) 4.0mm
Reusable 162-510-04



ROUND LOOP ELECTRODE (7.0cm) 6.0mm
Reusable 162-510-06



ROUND LOOP ELECTRODE (7.0cm) 8.0mm
Reusable 162-510-08



ROUND LOOP ELECTRODE (7.0cm) 10.0mm
Reusable 162-510-10



PLATE ELECTRODE SMALL
Reusable 162-560



PLATE ELECTRODE LARGE
Reusable 162-570



LOOP ELECTRODES (4.0mm)

Raz Enterprises

P.O. BOX NO. 449.
SIALKOT-51310-Pakistan.
Phone: +92-52-3551857
Fax: +92-52-3558041
E-mail: Info@Raz-Enterprises.com



					
Size	1.0X0.4 cm	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm
Reusable	162-170-04	162-170-05	162-170-07	162-170-08	162-170-10





					
Size	1.0X0.5 cm	1.0X0.7 cm	1.0X0.8 cm	1.0X1.0 cm	1.0X1.5 cm
Reusable	162-175-05	162-175-07	162-175-08	162-175-10	162-175-15



						
Size	1.5X0.5 cm	1.5X0.7 cm	1.5X0.8 cm	1.5X1.0 cm	1.5X1.2 cm	1.5X1.5 cm
Reusable	162-180-05	162-180-07	162-180-08	162-180-10	162-180-12	162-180-15



						
Size	2.0X0.8 cm	2.0X1.0 cm	2.0X1.2 cm	2.0X1.3 cm	2.0X1.5 cm	2.0X2.0 cm
Reusable	162-190-08	162-190-10	162-190-12	162-190-13	162-190-15	162-190-20

HIGH FREQUENCY INNOVATION



LAPAROSCOPIC ELECTRODES

WITHOUT SUCTION

L-HOOK ELECTRODE (36cm)
150-200



J-HOOK ELECTRODE (36cm)
150-210



NEEDLE ELECTRODE (36cm)
150-220



SPATULA ELECTRODE (36cm)
150-230



BALL ELECTRODE (36cm)
150-240



WITH SUCTION

L-HOOK ELECTRODE (36cm)
150-205



J-HOOK ELECTRODE (36cm)
150-215



NEEDLE ELECTRODE (36cm)
150-225



SPATULA ELECTRODE (36cm)
150-235



BALL ELECTRODE (36cm)
150-245





ARTHROSCOPIC ELECTRODES



DAGGER (12 cm)
150-260



DAGGER (18 cm)
150-265



HOOK (12 cm)
150-270



HOOK (18 cm)
150-275



HOCKEY (12 cm)
150-280



HOCKEY (18 cm)
150-285



FRAZIER SUCTION TUBE 11 FRENCH (5 $\frac{1}{4}$ ", 13.5cm)
150-300



FRAZIER SUCTION TUBE 11 FRENCH (8", 20.5cm)
150-310



FRAZIER SUCTION TUBE 11 FRENCH (16", 40.6cm)
150-320





GRAVES VAGINAL SPECULUM
WITH SMOKE EVACUATION TUBE

LARGE	MEDIUM	SMALL
130-100	130-102	130-103



DISPOSABLE TUBE
LARGE 130-100-DT
MEDIUM 130-102-DT
SMALL 130-103-DT



DISPOSABLE TUBE
LARGE 130-117-DT
MEDIUM 130-118-DT



VU-MAX SPECULUM
WITH SMOKE EVACUATION TUBE

GRAVES (MEDIUM)	PEDERSON (MEDIUM)
130-105	130-115



DISPOSABLE TUBE
GRAVES (MEDIUM) 130-105-DT
PEDERSON (MEDIUM) 130-115-DT

VU-MORE (GRAVES) SPECULUM
WITH SMOKE EVACUATION TUBE

LARGE	MEDIUM
130-117	130-118



PEDERSON VAGINAL SPECULUM
WITH SMOKE EVACUATION TUBE

LARGE	MEDIUM	SMALL
130-200	130-202	130-203



DISPOSABLE TUBE
LARGE 130-200-DT
MEDIUM 130-202-DT
SMALL 130-203-DT

**ELECTROSURGICAL
INSTRUMENTS FOR
GYNECOLOGY**



ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY

Raz Enterprises

P.O. BOX NO. 449.
SIALKOT-51310-Pakistan.
Phone: +92-52-3551857
Fax: +92-52-3558041
E-mail: Info@Raz-Enterprises.com



COLLIN SPECULUM
WITH SMOKE EVACUATION TUBE
LARGE 130-111 **MEDIUM** 130-112 **SMALL** 130-113



CUSCO SPECULUM
WITH SMOKE EVACUATION TUBE
LARGE 130-108 **MEDIUM** 130-109 **SMALL** 130-110



WEISMAN SPECULUM
LEFT OPEN
LARGE 130-220 **MEDIUM** 130-224



WEISMAN SPECULUM
RIGHT OPEN
LARGE 130-221 **MEDIUM** 130-225

DRESSING FORCEPS
130-180 (21 cm)
130-181 (25 cm)



TISSUE FORCEPS
130-185 (21 cm)
130-186 (25 cm)



IRIS HOOK (25 cm)
130-175



EMMET HOOK (25 cm)
130-170



STRAIGHT HOOK (25 cm)
130-176



TWO-PRONG HOOK (25 cm)
130-178



THREE-PRONG HOOK (25 cm)
130-179





ELECTROSURGICAL INSTRUMENTS FOR GYNECOLOGY



LATERAL VAGINAL
RETRACTOR 130-120



KOGAN ENDOCERVICAL
SPECULUM



3mm 130-133



4mm 130-134



5mm 130-135

**Raz
Enterprises**

P.O. BOX NO. 449,
SIALKOT-51310-Pakistan.
Phone: +92-52-3551857
Fax: +92-52-3558041
E-mail: Info@Raz-Enterprises.com



LATERAL RETRACTOR
(WIDE OPEN) 130-122



RIGBY RETRACTOR
130-150




GELPI RETRACTOR
130-140

HIGH FREQUENCY INNOVATION



ELECTROSURGICAL INSTRUMENTS

FOR GYNECOLOGY



SCHROEDER TANACULUM (24cm)
130-160




SCHROEDER TANACULUM (24cm)
130-161



SINGLE TOOTH TANACULUM (24cm)
130-162



CAMPION FORCEPS (22cm)
130-168



RING FORCEPS (22cm)
130-165



LASER SURGERY INSTRUMENTS FOR GYNECOLOGY



GRAVES VAGINAL SPECULUM
WITH SMOKE EVACUATION TUBE

LARGE 140-100
MEDIUM 140-102
SMALL 140-103




PEDERSON VAGINAL SPECULUM
WITH SMOKE EVACUATION TUBE

LARGE 140-200
MEDIUM 140-202
SMALL 140-203




COLLIN SPECULUM
WITH SMOKE EVACUATION TUBE

LARGE 140-111
MEDIUM 140-112
SMALL 140-113



**LATERAL VAGINAL
RETRACTOR**
140-120



**KOGAN ENDOCERVICAL
SPECULUM**



3mm 140-133



4mm 140-134



5mm 140-135



LASER SURGERY INSTRUMENTS

FOR GYNECOLOGY



GELPI RETRACTOR
140-140



RIGBY RETRACTOR
140-150



SCHROEDER TENACULUM
140-160



IRIS HOOK
140-175



EMMETT HOOK
140-170



STRAIGHT HOOK
140-176



TWO-PRONG HOOK
140-178



THREE-PRONG HOOK
140-179